



# GEMINI

## Urodynamics Product Catalog

An expanded market for urology consumables and equipment puts patients and practices first.



# Contents

Atmos Air-Charged Catheters .....	page 4
Atmos Air-Charged Pressure Sensors .....	page 6
Urodynamics Systems .....	page 8
Uroflowmeters .....	page 12
Electrodes .....	page 13
Tubing & Extension Lines .....	page 14

# Designed for practitioners and their patients

## Atmos Air-Charged Urodynamics Catheters

Premium performance, easy administration, and elevated patient comfort

With the thinnest diameter of urodynamics catheters available, the Atmos line of air-charged catheters are designed for maximum patient comfort. Sterile, disposable, and color-coded for ease of use, Atmos gives you options for all your diagnostic needs.



Product Number	Description	Compatibility	Quantity
1A-AC	Single Sensor Abdominal Air-Charged Urodynamics Catheter	6FR, 25" long	10
1B-AC	Single Sensor Bladder Air-Charged Urodynamics Catheter	6FR, 25" long	10
1CB-AC	Single Sensor Coudé Bladder Air-Charged Urodynamics Catheter	6FR, 25" long	10
2BU-AC	Dual Sensor Bladder/Urethral Air-Charged Urodynamics Catheter	6FR, 25" long	10
1RB-AC	Single Sensor Radiopaque Bladder Air-Charged Urodynamics Catheter	6FR, 25" long	10
2RBU-AC	Dual Sensor Radiopaque Bladder/Urethral Air-Charged Urodynamics Catheter	6FR, 25" long	10



Proven Accuracy



Advanced Structure



Matched Performance



Smaller Diameter




## Single Sensor Abdominal Catheter

Conduct accurate abdominal pressure readings rectally or vaginally.




## Single Sensor Bladder Catheter

Provides accurate bladder pressure readings for standard urodynamic procedures.



## Single Sensor Coudé Bladder Catheter

With a Coudé tip, this model is designed to provide accurate bladder pressure readings for patients with partially blocked urethras or enlarged prostates.



## Dual Sensor Bladder / Urethral Catheter

Two pressure channels allow you to measure both bladder and urethral pressures simultaneously. This model grants physicians an effective way of monitoring patients Urethral Pressure Profiles (UPP) and diagnosing urethral instability.

# Atmos Air-Charged Pressure Sensors

Reusable pressure sensors connect Atmos Air-Charged Catheters to your urodynamics system. Color fittings coordinate with connector on catheter for fast and accurate set-up.



Product Number	Description	Compatibility	Quantity
PS1011-V2M5	Atmos Pressure Sensor - Vesical 2.5m	MMT Systems except Harmony	1
PS1011-A2M5	Atmos Pressure Sensor - Abdominal 2.5m	MMT Systems except Harmony	1
PS1011-U2M5	Atmos Pressure Sensor - Urethral 2.5m	MMT Systems except Harmony	1
PS1021-V0M5	Atmos Pressure Sensor - Vesical 0.5m	MMT Harmony	1
PS1021-A0M5	Atmos Pressure Sensor - Abdominal 0.5m	MMT Harmony	1
PS1021-U0M5	Atmos Pressure Sensor - Urethral 0.5m	MMT Harmony	1
PS9000	MMT Holder for Atmos Pressure Sensors	Holds 4 Pressure Sensors	1

Product Number	Description	Compatibility	Quantity
PS1131-V2M5	Atmos Pressure Sensor - Vesical 2.5m	Laborie Aquarius / Goby	1
PS1131-A2M5	Atmos Pressure Sensor - Abdominal 2.5m	Laborie Aquarius / Goby	1
PS1131-U2M5	Atmos Pressure Sensor - Urethral 2.5m	Laborie Aquarius / Goby	1
PS1131-V0M5	Atmos Pressure Sensor - Vesical 0.5m	Laborie Roam	1
PS1131-A0M5	Atmos Pressure Sensor - Abdominal 0.5m	Laborie Roam	1
PS1131-U0M5	Atmos Pressure Sensor - Urethral 0.5m	Laborie Roam	1
PS1221-V2M5	Atmos Pressure Sensor - Vesical 2.5m	MMS Solar System	1
PS1221-A2M5	Atmos Pressure Sensor - Abdominal 2.5m	MMS Solar System	1
PS1221-U2M5	Atmos Pressure Sensor - Urethral 2.5m	MMS Solar System	1
PS1221-V0M5	Atmos Pressure Sensor - Vesical 0.5m	MMS Wireless Module	1
PS1221-A0M5	Atmos Pressure Sensor - Abdominal 0.5m	MMS Wireless Module	1
PS1221-U0M5	Atmos Pressure Sensor - Urethral 0.5m	MMS Wireless Module	1

# Better urodynamics systems for your practice

## Symphony SC

Fully loaded system for complex urodynamics, UPP and video UDS

Symphony SC is the most advanced, innovative, complete urodynamics system on the market. This all-in-one platform performs video UDS and UPP. The system is easy to use and configure with dual touch screen interfaces and adjustable ergonomic workstation.

Advanced features include valsalva leak point pressure camera, integrated 3D bladder scanner and infusion weight.



## Symphony ML

Premium system designed for advanced urodynamics

Symphony ML is an advanced, state-of-the-art trolley system uniquely tailored to meet the needs of medium to high-volume practices. With a touch screen interface and adjustable ergonomic workstation, Symphony ML is upgradeable to perform UPP studies and offers additional advanced UDS options making it the perfect fit for practices looking for more functionality.





Intuitive & Modern



Smarter Reporting



User-Friendly Software



On-Demand Availability

## Symphony EL

Economical system for everyday urodynamics

Symphony EL is a compact, cost-effective trolley system designed for low to medium volume practices. With a built-in monitor and adjustable trolley, Symphony EL is ideal for practices looking to perform comprehensive urodynamic procedures.



## Melody

Simple, compact system designed for portable urodynamics

Melody is a complete, cost-effective urodynamics system mounted to an IV-pole for true portability. Wireless and modular, Melody uses Bluetooth technology to connect to a laptop to perform a full range of urodynamic procedures.





# Advanced features to elevate your practice



## Valsalva Leak Point Pressure Camera

Flexible camera allows easy visualization of urethral leaks in real-time with automated reporting.

Included with Symphony SC; Optional add-on with Symphony ML



## Integrated 3D Bladder Scanner

Portable, intuitive and accurate scanning provides real-time PVR information in seconds with automated reporting.

Included with Symphony SC; Optional with Symphony ML and Danflow Wave Flowmeter



## Infusion Weight Sensor

Precise calculated infusion volume based on innovative IV scale.

Included with Symphony SC; Optional with Symphony ML

## Compare Models

	Symphony SC	Symphony ML	Symphony EL	Melody
Complex CMG	✓	✓	✓	✓
Voiding Pressure Study	✓	✓	✓	✓
EMG	✓	✓	✓	✓
Wireless Uroflow	✓	✓	✓	✓
UPP Study	✓	Optional		
Manual UPP	✓	Optional	Optional	Optional
Video UDS	✓			
Video Leak Point Camera	✓	Optional		
Infusion Weight Sensor	✓			
Integrated Bladder Scanner	✓	Optional		
System Structure	Advanced Electronic Trolley	Advanced Electronic Trolley	Adjustable Trolley	Pole
Computer	Integrated Computer	Integrated Computer	Integrated Computer	Laptop
Monitor	2 Touchscreens	Touchscreen	Monitor	Laptop
Channels	16	8	5	5
Warranty	4 years	3 years	3 years	2 years
Procedure Table	Optional	Optional	Optional	Optional
DICOM EMR	Optional	Optional	Optional	Optional
HL7 EMR	Optional	Optional		

# Everything you need in one place

## Danflow Wave

### Intelligent, reliable uroflowmetry

Danflow Wave is a smart flowmeter that wirelessly connects to a laptop for quick, easy and accurate uroflowmetry.

- Recording starts automatically when your patient starts voiding
- One-hand height adjustable stand easily accommodates standing patient or placed under commode
- Waterproof weight cell improves hygiene and simplifies cleaning
- Seamless EMR, HL7 and DICOM integration
- Optional: Integrated 3D Bladder Scanner, Laptop



## Danflow Wave CU

### Compact, wireless uroflowmetry


Danflow Wave CU is a smart flowmeter for real-time, easy and accurate uroflowmetry at the touch of a button.

- Recording starts automatically when your patient starts voiding
- One-hand adjustable stand easily accommodates standing patient or placed under commode
- Waterproof weight cell improves hygiene and simplifies cleaning
- Built-in printer for real-time reports



## Electrodes

Our assortment of disposable Electrodes will meet the needs of your urological procedures with strong signal quality, comfortable skin adhesion for your patients, and attentive sanitation.



Product Number	Description	Quantity
CND1720003	Cleartrode Repositionable ECG Electrodes Snap Patch - 100 Sets	300
RA36-03	Three 36"/91.4cm snap lead wires - Black, Green, White	3
RA36-03II	Three 36"/91.4cm snap lead wires - Black, Red, White	3
CND1741003	Neotrode II ECG Electrode, Tape with Wire Attached, 3/Pouch - 100 Sets	300

## Pump tubing & Extension lines

Compatible with a variety of industry-leading brands and devices, our selection of pump tubes are designed to attain precise results while using the equipment and tools you're used to.



Product Number	Description	Compatibility	Quantity
UD-A36205	Pump Tubing 3,2x6,4mm	Laborie Legacy	25
MMT-1173	Urodynamic set for pump	MEDKONSULT	25
UD-A46215	Pump Tubing 4x6mm	Mediwatch	100
UD-A58245	Pump Tubing 5x8mm	MMS	25
UD-AELN	Transparent Extension Line (1500mm)	-	50



Gemini Medical Technologies LLC  
123 NW 13th Street, Suite 212  
Boca Raton, FL 33432

Phone: +1.833.943.6464  
Email: [info@geminimedtech.com](mailto:info@geminimedtech.com)  
[geminimedtech.com](http://geminimedtech.com)

LBL-MM-AC-4008 Rev2