

GEMINI

Smart bladder scanner made simple

Delivering precision and accuracy
in every scan



VITAcon.

Better bladder scanners for your practice

Vitascan cVue

Built for a streamlined user experience with unmatched accuracy

The Vitascan cVue is a simple, purpose-built system with advanced detection algorithm technology that gives you confidence in every scan.

- Automatically calculates bladder volume, provides bladder contours, and adjusts to avoid interference from the pubic bone
- Optional thermal printer mounted to the cart for easy access enables instant, real-time results
- Secure data transfer through Wi-Fi via SSID or USB
- 10" touchscreen display with adjustable tilt and height
- 5-wheel cart system for improved stability



✔ Intuitive and modern

✔ Easy to get started

✔ Reliably precise

✔ PVR measurements

Vitascan LT-C6

Designed for seamless integration and flexible data management

The Vitascan LT-C6 is a Windows®-powered system offering seamless EMR integration, remote support capabilities, and advanced data management options.

- Automatically calculates bladder volume and visualizes anomalies that may affect volume calculation
- Save studies in real-time locally or to a network in standard file formats for advanced EMR integration
- Automatically measures Post-Void Residual (PVR) with pubic bone detection technology
- 10" antimicrobial medical grade tablet
- 5-wheel cart system for improved stability



Choose the right bladder scanner that meets your practices' needs

FEATURES		CVUE	LT-C6
Real-time bladder image		Multi-plane scan from one position on the body (3D) with exact measurements of bladder volume and contours with ability to manually adjust as needed	
Bladder volume measurement		+/- 7.5% accuracy, includes Post-Void Residual (PVR)	
Screen		10.1" color touchscreen, supports latex gloves	9.7" touchscreen
Operating system		Vitacon's proprietary software	Windows® OS (10-11)
Printing		Through networked devices or optional on-cart thermal printer	Through networked devices
EMR integration		Connect to network wirelessly via SSID	Saves to network in real-time (PDF, JPEG, DICOM)
Barcode scanner		Optional	-
Volume range		0-1000 mL	
Probe frequency		3.45 MHZ	3.5 MHZ
Probe weight		0.84 lbs	0.8 lbs
Scanning method		Sector, 180 degrees	
Rotating positions		12 or 24 rotating positions	
Sweep angle		130 degrees	
Console dimensions		10.35 x 8.4 x 4.92 in	7.3 x 1.8 x 1.8 in
Warranty	Components	5 years	
	Computer	5 years	
	Accessories	5 years	
	Downtime	Advanced exchange	
	Remote service/support	3 years	
	Accidental drop coverage for probes	2 years	

Get more info at:
geminimedtech.com/bladder-scanners



Gemini Medical Technologies LLC

Phone: +1.833.943.6464

Email: info@geminimedtech.com

geminimedtech.com

123 NW 13th Street, Suite 212 Boca Raton, FL 33432