

GEMINI

Atmos™ air-charged urodynamics catheters



A better choice for physicians and patients



Discover Atmos - the first innovation in air-charged urodynamics catheters in twenty years. Designed for accurate results, efficient usability, and maximum comfort, Atmos matches the standards of today's leading catheter while enhancing the diagnostic experience for both patients and practitioners through innovative new features:

- Flexible design for a relaxed fit
- The smallest diameter on the market
- Compatible with your existing equipment
- Consistent and reliable results

Proven accuracy and performance

Tested with 151 catheters at 1,991 pressure points, Atmos delivered near perfect statistical agreement-even in equipment calibrated to a T-DOC® Air-Charged™ Urodynamics Catheter.²

STATISTICAL APPROACH	STATISTIC	95% CONFIDENCE INTERVAL	P
Lin's Concordance Correlation Coefficient Pc	.999	.999, >.999	<.001
Pearson's r	1.00		<.001
Bland & Altman's Limits of Agreement	-0.897	-2.884, 1.090	



Atmos

Everything you need in one place

We provide an assortment of quality essential urodynamic consumable items to simplify procurement creating efficiency for your practice.

Consumables

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	BOX QTY
1A-AC	Single Sensor Abdominal Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1B-AC	Single Sensor Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1CB-AC	Single Sensor Coudé Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
2BU-AC	Dual Sensor Bladder/Urethral Air-Charged Urodynamics Catheter	6 FR, 25" long	10

Pump Tubing and Extension Lines

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	BOX QTY
MMT-1173	Urodynamic Set for Pump	MEDKONSULT Medical Technology Equipment	25
UD-A36205	Pump Tubing, 3,2x6,4 mm	Laborie Legacy Equipment	25
UD-A58245	Pump Tubing, 5x8	Medical Measurement System Equipment	25
UD-AELN	Transparent Extension Line, 1500 mm	-	50

Electrodes

PRODUCT NUMBER	DESCRIPTION	BOX QTY
CND1720003 - 100	Cleartrode Repositionable ECG Electrodes Snap Patch (100 Sets)	300
CND1741003 - 100	Neotrode II Pediatric/Neonatal ECG Electrode, Tape with Wire Attached, 3/Pouch (100 Sets)	300

Choose Atmos

With access to premium products at a responsible price, you have the power to choose what's best for your program.

Visit geminimedtech.com to request a sample or place an order.

GEMINI | **Atmos**

Gemini Medical Technologies LLC
8211 Ulmerton Road
Largo, FL 33771

Phone: +1.833.943.6464
Email: info@geminimedtech.com
geminimedtech.com

* GEMINI and ATMOS are trademarks of Gemini Medical Technologies LLC. All Rights Reserved. Product names, brands and other trademarks featured or referred to within the GEMINI website and marketing documents are the property of their respective trademark owners

* T-DOC® is a registered trademark of Laborie Medical Technologies Corp

LBL-MM-AC-4007_Rev.4

