

GEMINI

Urodynamics Product Catalog

An expanded market for urology consumables and equipment puts patients and practices first



Your trusted partner in urodynamics

Gemini provides advanced urodynamic solutions and tailored service to improve best practices and quality care.

By understanding and fulfilling the unique needs of today's physicians and their patients, we believe we can build long-term partnerships that will shape the future of urologic diagnostics and generate industry-wide prosperity.




Table of Contents

Catheters	4
UDS Systems	10
Uroflowmeters	14
Bladder Scanners	16
Examination and Procedure Chairs	18



Atmos[®] Air-Charged Urodynamic Catheters



With the thinnest diameter of urodynamics catheters available, the Atmos line of air-charged catheters are designed for maximum patient comfort. Sterile, disposable, and color-coded for ease of use, Atmos gives you options for all your diagnostic needs.



Smallest diameter

At just 6FR outer diameter, the Atmos family of catheters have the smallest diameter of any urodynamics catheter available



Proven accuracy

Tests show that Atmos catheters deliver reliable and accurate results, matching standards set by ICS Guidelines



Flexible for comfort

Unlike other catheters that remain rigid, Atmos' medical-grade polyurethane softens inside the body and bends as patients move



Full range diagnostics


Measure and record the full spectrum of bladder, urethral, and abdominal pressures all with Atmos

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	QUANTITY
1A-AC	Single Sensor Abdominal Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1B-AC	Single Sensor Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
1CB-AC	Single Sensor Coudé Bladder Air-Charged Urodynamics Catheter	6 FR, 25" long	10
2BU-AC	Dual Sensor Bladder/Urethral Air-Charged Urodynamics Catheter	6 FR, 25" long	10

Single Sensor Abdominal Catheter

Conduct accurate abdominal pressure readings rectally or vaginally.

Product number: 1A-AC



Single Sensor Bladder Catheter

Provides accurate bladder pressure readings for standard urodynamic procedures.

Product number: 1B-AC



Single Sensor Coudé Bladder Catheter

With a Coudé tip, this model is designed to provide accurate bladder pressure readings for patients with partially blocked urethras or enlarged prostates.


Product number: 1CB-AC



Dual Sensor Bladder/Urethral Catheter

Two pressure channels allow you to measure both bladder and urethral pressures simultaneously. This model grants physicians an effective way of monitoring patients Urethral Pressure Profiles (UPP) and diagnosing urethral instability.

Product number: 2BU-AC



Atmos Air-Charged Pressure Sensors

Conduct accurate abdominal pressure readings rectally or vaginally.




Plug and play

The colors of the luer locks make attaching the catheters with the same color simple



Compatible

We have sensors compatible with every major UDS system in the US market



"Gemini offers a great product at a very competitive price. We are delighted with their service. When I need something shipped quickly, they are all too happy to help—just a great partner to do business with."

CENTRAL FLORIDA UROGYNECOLOGY

PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	QUANTITY
PS11-01V25	Atmos Pressure Sensor - Vesical 2.5m	MMT Systems except Harmony	1
PS11-01A25	Atmos Pressure Sensor - Abdominal 2.5m	MMT Systems except Harmony	1
PS11-01U25	Atmos Pressure Sensor - Urethral 2.5m	MMT Systems except Harmony	1
PS11-02V05	Atmos Pressure Sensor - Vesical 0.5m	MMT Harmony	1
PS11-02A05	Atmos Pressure Sensor - Abdominal 0.5m	MMT Harmony	1
PS11-02U05	Atmos Pressure Sensor - Urethral 0.5m	MMT Harmony	1
PS9000	MMT Holder for Atmos Pressure Sensors	Holds 4 Pressure Sensors	1
PS11-11V25	Atmos Pressure Sensor - Vesical 2.5m	Laborie Aquarius/ Goby	1
PS11-11A25	Atmos Pressure Sensor - Abdominal 2.5m	Laborie Aquarius/ Goby	1
PS11-11U25	Atmos Pressure Sensor - Urethral 2.5m	Laborie Aquarius/ Goby	1
PS11-11V05	Atmos Pressure Sensor - Vesical 0.5m	Laborie Roam	1
PS11-11A05	Atmos Pressure Sensor - Abdominal 0.5m	Laborie Roam	1
PS11-11U05	Atmos Pressure Sensor - Urethral 0.5m	Laborie Roam	1
PS11-21V25	Atmos Pressure Sensor - Vesical 2.5m	MMS Solar System	1
PS11-21A25	Atmos Pressure Sensor - Abdominal 2.5m	MMS Solar System	1
PS11-21U25	Atmos Pressure Sensor - Urethral 2.5m	MMS Solar System	1
PS11-21V05	Atmos Pressure Sensor - Vesical 0.5m	MMS Wireless Module	1
PS11-21A05	Atmos Pressure Sensor - Abdominal 0.5m	MMS Wireless Module	1
PS11-21U05	Atmos Pressure Sensor - Urethral 0.5m	MMS Wireless Module	1

Pump Tubing and Extension Lines

Compatible with a variety of industry-leading brands and devices, our selection of pump tubes are designed to attain precise results while using the equipment and tools you're used to.



PRODUCT NUMBER	DESCRIPTION	COMPATIBILITY	QUANTITY
UD-A36205	Pump Tubing 3,2x6,4mm	Laborie Legacy	25
MMT-1173	Urodynamic set for pump	MEDKONSULT	25
UD-A58245	Pump Tubing 5x8mm	MMS	25
UD-AELN	Transparent Extension Line (1500mm)	-	50

"We have had nothing but great results with Gemini! It is super easy to place orders and our rep is always happy to help and answer any questions we have! The affordable pricing and great customer service makes Gemini a no-brainer!"

NORTH AUSTIN UROLOGY SPECIALISTS

Electrodes

Compatible with a variety of industry-leading brands and devices, our selection of pump tubes are designed to attain precise results while using the equipment and tools you're used to.



PRODUCT NUMBER	DESCRIPTION	QUANTITY
CND1720003	Cleartrode Repositionable ECG Electrodes Snap Patch – 100 Sets	300
CND1741003	Neotrode II ECG Electrode, Tape with Wire Attached, 3/Pouch – 100 Sets	300

A one-stop-shop for urology practices

Gemini's lineup of urodynamic consumables and equipment are available and ready to ship to you. Is it time to switch up suppliers?





Urodynamic Systems



Focus on your patients, not fiddling with equipment. Our lineup of systems are designed to meet every practice's needs. All UDS systems are loaded with state-of-the-art software that pre-populates diagnosis and treatment options and are built to be reliable and highly compatible.



Intuitive and modern

Designed to meet every practice's needs with advanced features such as integrated post void bladder scan and leak point camera



Smarter reporting

State-of-the-art software with upgraded user experience using pre-populated diagnosis and treatment options



On-demand availability

Standardized system configurations with quick delivery turnarounds



User-friendly software

UDMvision's clean interface provides automated study phase navigation

MODEL	DESCRIPTION
Symphony SC	Fully loaded system for advanced UDS, Video UDS and UPP
Symphony ML	Premium system, upgradable to advanced UDS and UPP
Symphony EL	Economical system designed for robust UDS
Melody	Compact system designed for portable UDS

Symphony SC

The most advanced, comprehensive system on the market. Dual touchscreen interfaces and an adjustable ergonomic workstation make it easy and efficient to navigate complex urodynamics procedures.



Symphony ML

A premium trolley system designed to meet the needs of medium to high-volume practices seeking advanced functionality. Includes a touchscreen interface and an adjustable ergonomic workstation.



Symphony EL

A streamlined system designed for practices needing an efficient, cost-effective system for everyday use. With a built-in monitor and adjustable trolley, it is compact and easy to maneuver.



Melody

A straightforward, compact system designed for practices seeking a portable solution. Mounted on an IV-pole, it uses radiofrequency technology to connect seamlessly to the included pre-configured laptop.



Urodynamic Systems Add-Ons

Integrated 3D Bladder Scanner

Portable, intuitive and accurate scanning provides real-time PVR information in seconds with automated reporting.



Infusion Weight Sensor

Precise calculated infusion volume based on innovative IV scale.



FEATURES	SYMPHONY SC	SYMPHONY ML	SYMPHONY EL	MELODY
Complex CMG	✓	✓	✓	✓
Voiding Pressure Study	✓	✓	✓	✓
EMG	✓	✓	✓	✓
Wireless Uroflow	✓	✓	✓	✓
Urethral Pressure Profile Study	✓	Optional		
Manual UPP	✓	Optional	Optional	Optional
Video UDS	✓			
Video Leak Point	✓	Optional		
Infusion Weight	✓			
Integrated Bladder Scanner	✓	Optional		
System Structure	Advanced Electronic Trolley	Advanced Electronic Trolley	Adjustable Trolley	Pole
Computer	Integrated Computer	Integrated Computer	Integrated Computer	Laptop
Monitor	2 Touchscreens	Touchscreen	Monitor	Laptop
Channels	16	8	5	5
Warranty	4 years	3 years	3 years	2 years
Procedure Table	Optional	Optional	Optional	Optional
DICOM EMR	Optional	Optional	Optional	Optional
HL7 EMR	Optional	Optional		



Uroflowmeters



Our smart flowmeters start recording automatically when your patient starts voiding. Both models come with a one-hand height adjustable stand to accommodate a standing patient or being placed under a commode. They're also easy to clean and calibrate themselves when a beaker is placed.



Waterproof

Waterproof weight cell improves hygiene and simplifies cleaning



Automatic

Recording starts automatically when your patient starts voiding



Easily adjustable

One-hand adjustable stand easily accommodates standing patient or placed under commode



Self calibrates

Automatically calibrates when a beaker is placed, ensuring accurate measurements with various beaker types

MODEL	DESCRIPTION
Danflow Wave	Connects wirelessly to a laptop or tablet for quick, easy and accurate uroflowmetry
Danflow Wave CU	Connects wirelessly to a thermal printer for real-time, easy and accurate uroflowmetry
Danflow Wave with Laptop and Integrated Bladder Scanner	Connects wirelessly to a laptop or tablet for quick, easy and accurate uroflowmetry with PVR.

Danflow Wave CU

Ideal for fast, standalone uroflowmetry testing. No need for a laptop. The thermal printer provides immediate zebra printouts of the test for analysis.



Danflow Wave

Connects to your laptop wirelessly with reporting that includes display graph and nomogram. Can connect with your EMR to bring in patient information into reports.



Danflow Wave with Laptop and Integrated Bladder Scanner

Bundle a laptop and bladder scanner to include post-void residual (PVR) measurement along with uroflowmetry. Can combine PVR and uroflowmetry in the same report.



VITACON.

Bladder Scanners



Get better, faster visualization of the bladder with our easy-to-use bladder scanners. Both systems can identify, track, and visualize bladder anomalies as you scan with real-time images with a fully-automated 3D output of the bladder and its volume.

✓ **Intuitive and modern**
The touchscreen interface is simple with large button displays for simple scanning

✓ **Reliably precise**
Automatically calculates bladder volume with +/- 7.5% accuracy

✓ **Easy to get started**
Boot it up and start scanning with data wirelessly sent to your EMR

✓ **Measures PVR**
Assess urinary retention without the need for catheterization

MODEL	DESCRIPTION
Vitascan cVue	Built for a streamlined user experience with unmatched accuracy
Vitascan LT – C6	Designed for seamless integration and flexible data management

Vitascan cVue

This scanner's advanced detection algorithm automatically calculates bladder volume, provides bladder contours, and adjusts to avoid interference from the pubic bone.



Vitascan LT-C6

This Windows®-powered system provides more opportunities for EMR integration, remote support, and advanced data management.

Examination and Procedure Chairs



Our modern and versatile chairs are designed specifically to meet the needs of urodynamics, urology, urogynecology, and OB/GYN practices. Their state-of-the-art technology and exceptional design is built to last and work for many procedures.

✓ Easily adjustable

Programmed with memory positions that can be easily switched by hand or foot controls

✓ Easy to use

Large seat cut-out makes access easy and fits with optional uroflow funnel holder

✓ Ergonomic

Comfortable seating position and two leg support options both include handles

✓ Strong and safe

Thick upholstery and stable base that supports at least 550 lbs

MODEL	DESCRIPTION
arco-matic® 300M	Adjusts to horizontal, lithotomy, high lithotomy, and seated positions automatically
medi-matic®	Meets ADA surgical standards and provides further patient accessibility options

arco-matic® 300 M

This chair is highly versatile for a range of examinations and procedures. It is available in six colors and with the choice of foot supports or göpel type leg supports.



medi-matic®

This chair can be used for outpatient surgeries such as transurethral resection of the prostate and vasectomy, in addition to exams and procedures. It supports more weight, has a lower starting height, and more memory positions.



Gemini Medical Technologies LLC

Phone: +1.833.943.6464

Email: info@geminimedtech.com
geminimedtech.com

8211 Ulmerton Road, Largo, FL 33771